## Multi-Lumen Balloon Applicator For Brachytherapy

### ENGLISH

**REF** B01245, B01256, B11245 & B11256

#### INSTRUCTIONS FOR USE

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

#### DESCRIPTION

The CONTURA® MLB Applicator consists of a multi-lumen catheter attached to an inflatable spherical balloon (Figure 1). Lumens are provided for attachment to commercially available HDR (High Dose Rate) remote afterloader equipment for passage of the radiation treatment delivery wire. Five treatment lumens are provided; one central lumen located along the long axis of the Applicator and four curved lumens symmetrically offset from the central lumen. A removable stiffening stylet is positioned in the central treatment lumen, and treatment lumens #1-4 have lumen caps. Two proximal ports are also provided with Luer-type connectors for balloon inflation/deflation and for application of intracavitary vacuum.

The CONTURA® MLB accessories provided for introduction and deployment include: trocar with split sheath, drainage catheter, three, 30 ml and one, 10 ml inflation syringes, #11 scalpel, contrast media tray, radiation lumen caps and labels (Figure 2).

#### Afterloader compatibility:

Varian	Nucletron
<ul style="list-style-type: none"> <li>VariSource ID™</li> <li>VariSource™ 200 Series</li> <li>GammaMedplus™ Series</li> <li>VariSource™ IX Series</li> </ul>	<ul style="list-style-type: none"> <li>microSelectron™ V2</li> <li>microSelectron™ Digital</li> <li>Flexitron™</li> </ul>

**Warning: The safety and effectiveness of the CONTURA® Applicator as a replacement for whole breast irradiation in the treatment of breast cancer has not been established.**

#### INDICATIONS FOR USE

The CONTURA® MLB Applicator is intended to provide brachytherapy when the physician chooses to deliver intracavitary radiation to the surgical margins following lumpectomy for breast cancer.

#### CONTRAINDICATIONS

- The Applicator is not intended for use in cavities that are too small, too large and/or of shapes unable to conform to the inflated balloon.
- The Applicator is not intended for use in patients with unusual anatomy including a highly curved rib structure and/or unequal amounts of tissue surrounding the cavity that may cause the CONTURA® balloon to be asymmetrical.

#### WARNINGS

- Use caution when positioning the trocar tip near the chest wall or skin margin to avoid unintended tissue damage.
- Do not fill the Applicator with more than 58 ml (B01245, B11245) or 108 ml (B01256, B11256) of fluid as overfilling may result in balloon rupture and/or device failure.
- The Applicator must be pre-tested before implantation. Do not use the balloon if it is not approximately symmetrical and/or any leakage is detected.
- The breast cavity must be imaged before implantation to insure the Applicator will fit appropriately. Do not use if the cavity is too small. Do not use if the balloon surface to skin surface distance is less than 5 mm unless the maximum skin dose is  $\leq 145\%$  of the prescription dose. Use of the appropriate off-set lumen(s) should be used to minimize the exposure to the skin.
- To insure appropriate treatment dose distribution, the Applicator must be imaged prior to delivering each fraction of radiation to confirm correct position, balloon volume, skin spacing and conformance.
- If excessive resistance is encountered when attempting to remove the Applicator from the patient, surgical removal is recommended.
- Contrast media concentrations of less than 10% are recommended to prevent dose attenuation.

- Non-ionic contrast media is recommended for patients who are allergic to iodine-based agents.
- This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminate period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.
- Do not resterilize. After resterilization, the sterility of the product is not guaranteed because of an indeterminate degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.

#### PRECAUTIONS

- The Applicator must be used only by physicians trained in catheter implantation, radiation treatment planning and delivery.
- Metal vascular and marking clips should not be used during the lumpectomy procedure to prevent potential abrasion or puncture of the CONTURA® balloon. Care should also be taken to direct suture knots and tails away from the cavity and whenever possible position tissue between the potential balloon surface and the tails.
- Store the CONTURA® Applicator at room temperature (25°C maximum).
- Care must be taken when handling and manipulating the CONTURA® balloon to prevent damage and foreign material contamination of the balloon surface.
- A scalpel should be used to incise the skin prior to inserting the trocar tip.
- Do not inject fluids into the Vacuum Port.
- Replace Luer caps and radiation lumen caps after use.
- Only clinical personnel trained in the operation of HDR afterloaders should deliver radiation using the Applicator.
- Verify that the appropriate afterloader connectors are available and function with the Applicator prior to treatment.
- Important: the Applicator must remain as straight as possible and free of sharp bends and kinks when connected to the HDR afterloader. Treatment length measurements should be obtained prior to each treatment fraction with the device in the same orientation as treatment. Orientation should remain consistent for each treatment fraction.
- Inspect package before use. Discard if seal is compromised or packaging is damaged.

#### COMPLICATIONS

Complications that may be associated with the use of the CONTURA® MLB Applicator are the same as those associated with the use of similar devices. These may include: erythema, catheter site drainage, breast pain, ecchymosis, breast fibrosis, telangiectasia, breast induration, breast seroma, breast edema, dry desquamation, dry skin, skin discoloration, parasthesia, axillary pain, fatigue, pruritis, breast retraction, nausea, skin irritation, moist desquamation, hematoma, rash, asymptomatic fat necrosis, breast infection, breast blister and lymphedema.

#### HOW SUPPLIED

The CONTURA® MLB Applicator and accessories are provided sterile and are intended for single patient use only.

#### DIRECTIONS FOR USE

#### PLACEMENT - Refer to Figures 1 & 2



Figure 1: CONTURA® APPLICATOR


**Figure 2: ACCESSORIES**

1. Use ultrasound to identify the lumpectomy cavity.
2. Open the **CONTURA®** MLB Applicator sterile package and remove the Applicator (A) and one 30 ml Syringe (B). Remove the Inflation Port Luer Cap (C) and inject 58 ml (B01245, B11245) or 108 ml (B01256, B11256) of sterile saline into the Applicator and inspect for leaks and symmetry. Discard Applicator if defective. Holding the Applicator by the connectors, with the balloon hanging vertically, completely withdraw the saline from balloon.
3. Prepare a maximum 3% contrast media/sterile saline solution in the Tray (D) provided.  
**Note:** Higher concentrations of contrast media may inhibit lumen marker visualization.
4. Determine the desired point on the breast surface for the insertion of the Applicator. Inject appropriate anesthetic to the skin and pathway to the lumpectomy cavity. Make a skin incision with the scalpel at the insertion point of sufficient length to fully insert the Trocar (F) tip. Dilate the skin incision, if desired. Advance the Trocar with Split Sheath (G) into the cavity. Remove the Trocar.
5. Attach a 30 ml syringe to the Drainage Catheter (H) and drain any fluids within the cavity by inserting the Drainage Catheter through the Split Sheath and suctioning. Remove the Drainage Catheter.
6. Insert the Applicator through the Split Sheath into the cavity. Remove the Sheath.
7. Align the Treatment Lumen #1 Stripe (I) on the catheter shaft with the skin incision.
8. Remove the stiffening Stylet (J) from the Central Treatment Lumen (K). Attach a red Treatment Lumen Cap (L).
9. Using the syringes provided, purge any air from the applicator and inflate the Applicator balloon with the contrast media solution to the desired fill volume. Purge any air from the fill syringes before attaching them to the Applicator.

Model	Desired balloon diameter	Approximate balloon fill volume
B01245	4 cm Dia. x 4.4 cm length	33 ml
B11245	5 cm Dia. x 4.4 cm length	58 ml
B01256	4.5 cm Dia. x 5.4 cm length	51 ml
B11256	6 cm Dia. x 5.4 cm length	108 ml

10. Replace the Luer Cap on to the Inflation Port (M).
11. Use ultrasound to confirm appropriate placement, volume and cavity conformance. Fluid and air surrounding the Applicator balloon may be aspirated with a 30 ml Syringe attached to the white Vacuum Port (N). The volume of the balloon may be adjusted through the blue Inflation Port (M). Replace Luer Caps when finished.
12. Confirm that the Treatment Lumen #1 Stripe is aligned with the skin incision.
13. Position catheter to minimize bending and to align marker with skin incision.
14. Apply a surgical dressing to the exit site with the catheter positioned to minimize bending.
15. Record the final balloon fill volume on the labels provided and attach to the patient's chart.

## **RADIATION DELIVERY - Refer to Figures 1 & 3**

**Figure 3: TREATMENT LUMEN ORIENTATION**

1. Important: the Applicator must remain as straight as possible and free of sharp bends and kinks when connected to the HDR afterloader. Treatment length measurements should be obtained prior to each treatment fraction with the device in the same orientation as treatment. Orientation should remain consistent for each treatment fraction.
2. CT imaging should be used in conjunction with commercially available treatment planning software to determine the appropriate treatment lumens, treatment dwell positions and dwell times for optimized radiation delivery of a prescribed dose to the targeted treatment volume.
3. Using Treatment Lumen #1 Stripe as a reference, note the orientation of the **CONTURA®** Applicator. Verify correct Applicator orientation, balloon position, balloon volume, skin spacing and conformance using imaging prior to delivery of each fraction of radiation. Adjust if necessary.
4. The Treatment Lumens are numbered '1', '2', '3', '4' and '5' and positioned as shown in Figure 3. Lumen number '1' corresponds to the offset lumen closest and parallel to the Treatment Lumen #1 Stripe (I) along the catheter shaft. Treatment Lumen number '5' corresponds to the central lumen. Remove the red caps and use commercially available connectors to attach the treatment lumens to the afterloader.  
**Note:** When using the GammaMedplus™ afterloader, the treatment lumens of the applicator must first be trimmed to length using the appropriate GammaMedplus™ length cutting gauge.
5. After each treatment replace the red caps and secure catheter shaft.

## **REMOVAL**

1. Remove the **CONTURA®** MLB Applicator by first attaching a syringe to the blue Inflation Port and deflating the balloon.  
**Note:** If difficulty is encountered deflating balloon with syringe:  
1) Re-attach syringe and securely rotate clockwise to completely activate the valve. If the balloon, still does not deflate, then  
2) Cut the blue Inflation Port tubing. The saline/contrast contents of the balloon will now drain from the end of the cut tubing  
3) If the balloon still cannot be deflated, insert a syringe with needle directly into the balloon. Deflate the balloon by withdrawing contents using the syringe.
2. Rotate and withdraw (unscrew) the Applicator from the cavity.

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2. CT

Rx Only

**Attention, See Instructions For Use**

Attention, voir le mode d'emploi  
Achtung! Siehe Gebrauchsanweisung  
Attenzione, leggere le istruzioni per l'uso  
Atención, consulte las instrucciones de uso  
Let op, zie Gebruiksaanwijzing  
Atenção, Consulte as Instruções de Utilização  
Προσοχή, βλ. Οδηγίες χρήσης  
Bemærk, Se Brugervejledning  
Obs! Se bruksanvisningen  
Huomio! Lue käyttöohjeet  
OBS! Les bruksanvisningen  
Uwaga: Należy zapoznać się z instrukcją użytkowania  
Figyelem, lásd a használati utasítást  
Pozor, viz návod k použití  
Dikkat: Kullanım Talimatlarına Bakınız

**Contents**

Contenu  
Inhalt  
Contenuto  
Contenido  
Inhoud  
Conteúdo  
Περιεχόμενα  
Indhold  
Innehåll  
Sisältö  
Innhold  
Zawartość  
Tartalom  
Obsah  
İçindekiler

LOT

**Lot Number**

Numéro de lot  
Chargennummer  
Numero di lotto  
Número de lote  
Lotnummer  
Número do Lote  
Αριθμός παρτίδας  
Lotnummer  
Lotnummer  
Eränumero  
Partinummer  
Numer partii  
Sorozatszám  
Číslo šarže  
Lot Numarası  
Lot

**Use By**

Date de péremption  
Verwendbar bis  
Data di scadenza  
Fecha de caducidad  
Uiterste gebruiksdatum  
Válido até  
Ημερομηνία λήξης  
Anvendes inden  
Användes före  
Käytettävä ennen  
Bruk før  
Termin ważności  
Szavatosság lejár  
Datum expirace  
Son Kullanma Tarihi

**Manufacturer**

Fabricant  
Hersteller  
Fabbicante  
Fabricante  
Fabrikant  
Fabricante  
Κατασκευαστής  
Producent  
Tillverkare  
Valmistaja  
Fabrikant  
Produdent  
Gyártó  
Výrobce  
Üretici

**Do Not Resterilize**

Ne pas restériliser  
Nicht erneut sterilisieren  
Non risterilizzare  
No volver a esterilizar  
Niet opnieuw steriliseren  
Não reesterilizar  
Μην επαναποστερώνετε  
Må ikke resteriliseres  
Får ej resteriliseras  
Ei saa steriloida uudestaan  
Må ikke steriliseres igjen.  
Nie należy sterylizować ponownie  
Újrasterilizálni tilos!  
Neresterilizujte  
Tekrar Sterilize Etmeyin

REF

**Catalogue Number**

Numéro de catalogue  
Katalognummer  
Numero di catalogo  
Número de referencia  
Catalogusnummer  
Número de Catálogo  
Αριθμός καταλόγου  
Katalognummer  
Katalognummer  
Luettelonumero  
Katalognummer  
Numer katalogowy  
Katalógusszám  
Katalogové číslo  
Katalog Numarası

**Single Use**

À usage unique  
Zur einmaligen Verwendung  
Monouso  
De un solo uso  
Voor eenmalig gebruik  
Não Reutilizar  
Μίας χρήσης  
Til engangsbrug  
För engångsbruk  
Kertakäyttöinen  
Engangsbruk  
Do jednorazowego użytku  
Egyszer használatos  
Pouze pro jednorázové použití  
Tek Kullanılmıktır

**STERILE R**

**Sterilized Using Irradiation**

Stérilisé par irradiation  
Durch Bestrahlung sterilisiert  
Sterilizzato mediante irradiazione  
Esterilizado mediante irradiación  
Gesteriliseerd met behulp van straling  
Esterilizado por irradiação  
Αποστειρωμένο με χρήση ακτινοβολίας  
Steriliseret ved brug af bestråling  
Steriliserad med strålning  
Steriloitu säteilyttämällä  
Sterilisert ved stråling  
Produkt sterylizowany przez napromienianie  
Sugárzással sterilizált  
Sterilizovano ozářením  
Radyoaktif Işınla Sterilize Edilmiştir

**Do Not Use If Package Is Damaged**

Ne pas utiliser si l'emballage est endommagé  
Nicht verwenden, wenn die Verpackung beschädigt ist  
Non utilizzare se la confezione è danneggiata  
No usar si el envase está dañado  
Niet gebruiken wanneer de verpakking beschadigd is  
Não utilizar se a embalagem estiver danificada  
Μη χρησιμοποιείτε το προϊόν εάν η συσκευασία του έχει υποστεί ζημιά  
Må ikke bruges, hvis emballagen er beskadiget  
Använd inte produkten om förpackningen skadats  
Ei saa käyttää, jos pakkaus on vaurioitunut  
Skal ikke brukes hvis pakningen er skadet  
Nie stosować, jeśli opakowanie jest uszkodzone  
Ne használja fel, ha a csomagolás sérült  
Nepoužívejte, pokud je obal poškozen  
Paket Hasarlıysa Kullanmayınız

**Bard and Contura are trademarks and/or registered trademarks of C.R. Bard, Inc. or an affiliate**

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**Copyright ©2012 C. R. Bard, Inc. All rights reserved. Printed in the USA.**

Copyright ©2012 C. R. Bard, Inc. Tous droits réservés. Imprimé aux États-Unis.  
Copyright ©2012 C. R. Bard, Inc. Alle Rechte vorbehalten. Gedruckt in USA.  
Copyright ©2012 C. R. Bard, Inc. Tutti i diritti riservati. Stampato negli Stati Uniti.  
Copyright ©2012, C. R. Bard, Inc. Reservados todos los derechos. Impreso en EE.UU.  
Copyright ©2012 C. R. Bard, Inc. Alle rechten voorbehouden. Gedrukt in de VS.  
Copyright ©2012, C. R. Bard, Inc. Todos os direitos reservados. Impresso nos E.U.A.  
Πνευματικά Δικαιώματα ©2012 C. R. Bard, Inc. Με την επιφύλαξη παντός δικαιώματος. Εκτυπώθηκε στις Η.Π.Α.  
Copyright ©2012 C. R. Bard, Inc. Alle rettigheder forbeholdes. Trykt i USA.  
Copyright ©2012 C. R. Bard, Inc. Alla rättigheter förbehållna. Tryckt i USA.  
Copyright ©2012 C. R. Bard, Inc. Kaikki oikeudet pidätetään. Painettu Yhdysvalloissa.  
Copyright ©2012 C. R. Bard, Inc. Alle rettigheter forbeholdt. Trykt i USA.  
Copyright ©2012 C. R. Bard, Inc. Wszystkie prawa zastrzeżone. Wydrukowano w USA.  
Copyright ©2012, C. R. Bard, Inc. Minden jog fenntartva. Nyomtatva az USA-ban.  
Copyright ©2012 C. R. Bard, Inc. Všechna práva vyhrazena. Vytisknuto v USA.  
Yapıt Hakkı ©2012 C. R. Bard, Inc. Tüm hakları saklıdır. ABD’de basılmıştır.  
© 2012 C. R. Bard, Inc.

 **Manufacturer:**  
**SenoRx, Inc.**  
1625 West 3rd Street  
Tempe, AZ 85281  
USA

TEL: 1-480-894-9515  
1-800-321-4254  
FAX: 1-480-966-7062  
1-800-440-5376  
[www.bardbiopsy.com](http://www.bardbiopsy.com)

**Authorized Representative in  
the European Community**  
Bard Limited  
Forest House  
Tilgate Forest Business Park  
Brighton Road, Crawley  
West Sussex  
RH11 9BP, UK

**BARID**  

---

**BIOPSY SYSTEMS**